

# RESCUE

## RESCUE RANDY MANIKINS

Developed for lifelike adult or juvenile patient handling, transportation, and extrication training. These manikins can be safely used in situations too hazardous or uncomfortable for human volunteers.

Made of tough cast vinyl with zinc-plated steel frame for longevity and rugged use. Both manikins are fully articulated for positioning in widely diversified training environments. The manikins can be seated in an automobile, positioned in a smoke-filled room for fire rescue training, and immobilized on a stretcher for mountain rescue instruction. They are also useful for simulating the removal of victims from utility poles, conduits, mine cave-ins, and tunnels, as well as in limitless other training procedures.

The represents 1,65 manikins tall and available unweighted at 23 kgs. or weighted to 48 kgs. with a new mixture of vinyl beads and foam to duplicate the weight and resiliency of a human body. Shipped in sweat pants.

### MSA-50

#### Unweighted "Rescue Randy".

Dimensions: 127 x 54 x 31 cm.

Sh. wt. 29.02 kgs.

### MSA-60

#### Weighted "Rescue Randy".

Dimensions: 127 x 54 x 31 cm.

Sh. wt. 56.69 kgs.

### MSA-70

#### Weighted "Rescue Randy".

Dimensions: 127 x 54 x 31 cm.

Sh. wt. 74.25 kgs.

### MSA-80

#### Weighted "Rescue Randy".

Sh. wt. 83.25 kgs.

Dimensions: 127 x 54 x 31 cm.

Available supplies

- Carry Storage Bag 0421336



MSA-70

## WATER RESCUE MANIKINS

These adult and adolescent models are made of durable plastic with rust-resistant skeletons, articulated joints, and optional CPR capabilities. As the manikins fill with water, they will sink up to the neck. To take them to the bottom, add five to ten pounds of weight (not supplied). A brick will do.

### MRA-100

#### Adult Water Rescue Manikin

Dimensions: 165 x 59 x 23 cm.

Sh. wt. 23.4 kgs.

### MRA-110

#### Adult CPR Water Rescue Manikin

Dimensions: 165 x 59 x 23 cm.

Sh. wt. 20.25 kgs.

### MRA-80

#### Adolescent CPR Water Rescue Manikin

Dimensions: 125 x 41 x 23 cm.

Sh. wt. 10 kgs.

### MRA-90

#### Adolescent Water Rescue Manikin

Sh. wt. 10 kgs.

Dimensions: 125 x 41 x 23 cm.



RESCUE ADULT



RESCUE ADOLESCENT