


NURSING CARE




SPC-01 Chloe Patient Care Simulator

A fully articulating patient care simulator with these airway management features: soft, lifelike airway with realistic tongue, trachea, and epiglottis; BVM with realistic chest rise; practice oral/nasal intubation; and pass NP/OP tubes. General care features include: stylish wig for hair care exercises and surgical draping; removable dentures; eyes in eye sockets for ophthalmic exercises; ear canal for otic drops and irrigation; injection arm for IV, IM, Sub Q, and intradermal exercises; pass NG and OG tubes; practice colostomy, ileostomy, and suprapubic care/irrigation exercises; interchangeable male/female catheterization; enema administration; decubitus ulcers; and ulcerated foot. Sh. wt. 36 kgs.




Breast Examination

Soft vinyl bra with 7 interchangeable breasts for palpation and examination. These breasts demonstrate chronic mastitis, benign growth, carcinoma and the "orange skin" effect, giant sarcoma, scirrhous carcinoma, the lymphatic drainage system, and movement of the mammary gland on the surface of the pectoralis major muscle.



Heart & Lung Sounds

Sensor network hidden beneath the skin. Hear the appropriate heart or lung sound as bell of stethoscope is moved across the front and back of the torso. Includes our Virtual Stethoscope™ with multiple heart and lung sounds. An external speaker plugs into the Virtual Stethoscope™ so a classroom can hear what the student hears.



Gynecologic Examination

Vaginal douching and PAP smear exercises. Bi-manual pelvic examination. Palpation of normal and pregnant uteri. Vaginal examination, including insertion of speculum. Visual recognition of normal and abnormal cervixes. Uterine sounding.



Blood Pressure Training arm

Adjustable systolic and diastolic pressures. Adjustable auscultation gap. Adjustable pulse rate. Display tracks cuff pressure.



Ulcerate foot



Interchangeable male/female catheterization



Injection arm for IV, IM sub q and intradermal

NURSING CARE

MEG-70

Complete Dual-Sex Geriatric Manikin

The Nursing Skills Manikin is ideal for preparing students for working with real-life geriatric patients, since it features an elderly appearance with skin wrinkles and folds. No other manikin on the market is so realistic and true-to-life. When compared to similar manikins, it features the best and most lifelike range of motion; Nursing moves like a person, allows realistic patient positioning, and features nonpinching joints. The overall appearance of the manikin is female, however, removal of the wig and attachment of the male genitals converts it for male catheterization and prostate examination simulations. Excellent for all levels of health care education... from nursing assistants to medical students. Ideal for all OBRA required training.

Dimensions: 97 x 51 x 28 cm.

Sh. 14.51 kgs.

Visual Inspection Features


- Decubitus Sacral Ulcer
- Cancerous Mole Comparison
- Dilated Pupil Comparison
- Elderly Appearance
- Skin Wrinkles and Folds

Fully Functional Patient Care Simulation

- Urinary Catheterization- Female and Male (Uncircumcised)
- Prostate Exam - Stage B - Male
- Ostomy Care - Ileostomy and Colostomy Lavage and Suctioning
- Tracheostomy Care - Lavage and Suctioning
- Gastrostomy Care, Lavage, Gavage
- **Optional** IV Arm Attachment
- Denture Removal - Upper and Lower
- Eye Irrigation
- Bed Baths
- Hair Care - Washing, Combing
- Clothing Changes
- Patient Transfer Techniques
- Enema Administration - Female
- Pap Smears and Douching
- Placement of Nasogastric Tubes
- Oral and basal Lavage, Gavage and Suctioning
- Intramuscular Injection Sites - Arm, Thigh, Buttock
- **Optional** Blood Pressure Arm Attachment
- Oral Hygiene
- Ear Irrigation and Hearing Aid Placement
- Finger and Toe Bandaging - Flexible and Individually Molded
- Bandaging and Wound Dressing
- Patient Positioning

Realistic Patient Transfer Procedures

- Trunk - Rotation, Hyperextension
- Shoulder - Abduction, Adduction, Rotation, Hyperextension
- Elbow - Extension, Flexion, Pronation, Supination
- Wrist - Flexion, Hyperextension, Radial Flexion, Ulnar Flexion
- Fingers - Abduction, Adduction, Flexion
- Neck - Rotation, Hyperextension, Lateral Flexion
- Hip - Abduction, Adduction, Rotation, Hyperextension
- Knee - Extension, Flexion
- Ankle - Eversion, Inversion, Dorsiflexion, Plantarflexion
- Toes - Abduction, Adduction, Flexion



MEG-70

MEG-90

Advanced Dual-Sex Geriatric Manikin

Same as model MEG-70, but also includes an IV Training Arm and a Blood Pressure Training Arm.

Dimensions: 97 x 51 x 28 cm.

Sh. wt. 19.59 kgs.

MEG-50

Basic Dual-Sex Geriatric Manikin

The basic Geri Manikin provides all the regular features of the Geri Complete Manikin, except internal fluid reservoirs have been eliminated for greater economy. Urinary catheterization cannot be performed on the Basic Geri.

Dimensions: 97 x 51 x 28 cm.

Sh. wt. 13.61 kgs.

Available Supplies

- Foley Catheter (16 FR. 5cc.) Pkg. of 1 0421000
- Foley Catheter (16 FR. 5cc.) Pkg. of 10 0421001

Accessories

BII-10

IV Arm

Sh. wt. 3.23 kgs.

SPB-40

Blood Pressure Arm

Sh. wt. 4.99 kgs.

AMS-100 (110 Volt) AMS-100E (220 Volt)

Speaker System for Blood Pressure Arm

- Soft Manikin Case 0421126
- Hard Manikin Case 0421124



MEG-90

NURSING CARE

MEK-60

Complete Dual-Sex Nursing Skillis Manikin

The Nursing Skills Manikin has all the features listed below without being age-specific. Molding and modeling techniques have given Nursing an unusually realistic appearance as a woman, or as a man when the wig is removed. Unique visual features reinforce the importance of inspecting patients and noting changes in their appearance. All joints provide maximum lifelike range of motion not found in any other manikin. Perfect manikin for all OBRA required training. Sh. wt. 15.28 kgs. Dimensions: 97 x 51 x 28 cm.

Visual Inspection Features include:

- Decubitus Sacral Ulcer - Stage 1
- Dilated Pupil Comparison
- Cancerous Mole Comparison
- Reddened Skin Folds

Complete Range of Motion include:

- TRUNK - Rotation, Hyperextension
- NECK - Rotation, Hyperextension, Lateral Flexion
- SHOULDER - Abduction, Adduction, Rotation, Hyperextension
- HIP - Abduction, Adduction, Rotation, Hyperextension
- ELBOW - Extension, Flexion, Pronation, Supination
- KNEE - Extension, Flexion
- WRIST - Flexion, Hyperextension, Radial Flexion, Ulnar Flexion
- ANKLE - Eversion, Inversion, Dorsiflexion, Plantarflexion
- FINGERS - Abduction, Adduction, Flexion (soft, lifelike material)
- TOES - Abduction, Adduction, Flexion (soft, lifelike material)

Lifeform Patient Care Simulations includes:

- Urinary Catheterization - Female and Male (Uncircumcised)
- Enema Administration - Female
- Prostate Exam, Stage B - Male
- Pap Smears and Douching
- Ostomy Care - Ileostomy and Colostomy Lavage and Suctioning
- Placement of Nasogastric Tubes
- Oral and Nasal Lavage, Gavage, and Suctioning
- Tracheostomy Care - Lavage and Suctioning
- Intramuscular Injection Sites - Arms, Thigh, Buttock
- Denture Removal - Upper and Lower
- Oral Hygiene
- Eye and Ear Irrigation
- Bed Baths
- Finger & Toe Bandaging - Flexible and Individually Molded
- Hair Care - Washing, Combing
- Bandaging and Wound Dressing
- Clothing Changes
- Patient Positioning and Transfer Techniques
- Hearing Aid Removal and Placement Techniques
- Optional IV Arm Attachment
- Optional Blood Pressure Arm

MEK-80

Advanced Dual-Sex Nursing Skillis Manikin

Same as model MEK-60, but also includes an IV Training Arm and a Blood Pressure Training Arm.

Dimensions: 97 x 51 x 28 cm. Sh. wt. 19.79 kgs.

MEK-40

Basic Dual-Sex Nursing Skillis Manikin

The basic Manikin provides all the regular features of the Complete Manikin, except internal fluid reservoirs have been eliminated for greater economy. Urinary catheterization cannot be performed on the Basic.

Dimensions: 97 x 51 x 28 cm. Sh. wt. 13.61 kgs.

BP-95

Mike and Michelle 5-Year-Old Pediatric Care Simulator

A sophisticated pediatric simulator of a 5-year-old for training in standard and advanced clinical procedures. The realistic face, body, and fully jointed elbows, wrists, knees, and ankles permit the student and teacher to simulate authentic nursing care and medical treatment of the young child. Features include: soft, lifelike facial skin; self-molded hair; eyes open and close in realistic eye sockets for ophthalmic procedures; fully articulated head and jaw with teeth and tongue; NG and otic exercises; bends 30° at waist; soft upper body skin over hard upper body for realistic feel; intramuscular injection sites in left deltoid and right upper thigh; TB test site; interchangeable male organ; tracheotomy opening; male and female catheterization; removable internal tanks; enema administration; and neck brace. Detachable at waist for easy storage. Comes with nylon carrying bag. Sh. wt. 10.43 kgs.



MEK-60

Available Supplies

- Foley Catheter (16 FR. 5cc.) Pkg. of 1 0421000

- Foley Catheter (16 FR. 5cc.) Pkg. of 10 0421001

Accessories

BII-10

IV Arm

Sh. wt. 3.23 kgs.

SPB-40

Blood Pressure Arm

Sh. wt. 4.99 kgs.

AMS-100 (110 Volt) AMS-100E (220 Volt)

Speaker System for Blood Pressure Arm

- Soft Manikin Case 0421126

- Hard Manikin Case 0421124

NURSING CARE



BI-68 Intravenous Training Arm Option (model BP-95)

Advanced injection training arm and hand is attachable to manikin. Sh. wt. 2.47 kgs.



BP-40 Nurse Training Baby

This advanced and incredibly realistic model features internal organs (heart, lung, intestine system, stomach, and bladder), a removable belly covering, interchangeable genital organs, soft ears, and fully movable head, arms, and legs. It allows students to practice both basic and medical infant care such as tracheotomy care, male and female catheterization, insertion of stomach tube, enema insertion, as well as femoral and gluteal intramuscular and subcutaneous injections. The fontanels can be palpated. Additionally, this model allows practice in dressing; changing; feeding simulation; navel care; mouth, ear, and nose cleaning; insertion of a rectal thermometer; and proper carrying techniques. Baby is 52 cm long and weighs 3.00 kgs. Comes with a carrying bag. Sh. wt. 4.50 kgs.



BP-62 Newborn Advanced Care Simulator

This newborn simulator looks just like a 0-8-week-old baby. The chubby face and body with beautifully detailed tiny, soft hands and feet make you want to hold and care for it. Features include: training arm and hand for intravenous access (also features radial and brachial arterial sites); external stoma sites with internal tanks; umbilical catheterization site with internal tank; temporal venous access sites with internal tank; soft and flexible facial skin; self-molded hair; realistic eyes; NG exercises to demonstrate tube feeding and gastric suction; simulated ear canal; tongue may be moved from side to side; soft arms and legs rotate within the torso body for lifelike feel and position; soft hands, feet, fingers, and toes; heel stick and finger prick technique; soft upper body skin over torso for "babylike" feel; bathing and bandaging activity; intramuscular injection site in upper thigh; TB test site; interchangeable male organ; urethral passage and bladder; male and female catheterization; removable internal tanks; enema administration; T-shirt and diapers; soft nylon carrying bag; and instruction manual.

Dimensions: 81 x 20 x 20 cm.

Sh. wt. 3.17 kgs.



BP-60 Newborn Care Simulator (not illustrated).

This newborn simulator looks just like a 0-8-week-old baby. The chubby face and body with beautifully detailed tiny, soft hands and feet make you want to hold and care for it. Features include: soft and flexible facial skin; self-molded hair; realistic eyes; NG exercises to demonstrate tube feeding and gastric suction; simulated ear canal; tongue may be moved from side to side; soft arms and legs rotate within the torso body for lifelike feel and position; soft hands, feet, fingers, and toes; heel stick and finger prick technique; soft upper body skin over torso for "babylike" feel; bathing and bandaging activity; intramuscular injection site in upper thigh; TB test site; interchangeable male organ; urethral passage and bladder; male and female catheterization; removable internal tanks; enema administration; T-shirt and diapers; soft nylon carrying bag; and instruction manual.

Dimensions: 81 x 20 x 20 cm.

Sh. wt. 3.00 kgs.